

Dog Ends™

Dog Ends™ assist in the rapid healing of tail tip injuries in dogs

- Patent pending design
- Ultralight
- Easy to apply

Dog Ends™ protect tail tip injuries from trauma and aid rapid healing. **Dog Ends™** fit over and overlap the end of the tail and are held in place by a small wrap of medical tape and the premium grade plastic mesh material takes the impact when the dog wags its tail. Unlike conventional fabric or solid dressings, the mesh will not absorb wound fluid and the open weave design ensures plenty of ventilation around the wound surface to promote clean healing. **Dog Ends™** have been fully tested to ensure the ideal combination of stretch and strength and are easy to apply using the specially designed Applicator. There is no need for bandages, although when necessary, **Dog Ends™** can be used to hold a simple wound dressing in place for first aid on a fresh wound. They can also be used preventatively to stop new wounds occurring.

Dog Ends™ come in a Starter Pack containing an Applicator, 10 mesh **Dog Ends™** and a roll of medical tape. (Refill packs of mesh, tape and Applicator are also available).



Hound Surrounds™

Hound Surrounds™ are a completely new design of Elizabethan collars for dogs

- Patent pending design
- Light and comfortable for the dog
- Flexible and strong

Hound Surrounds™ are primarily for bigger dogs, as it is the larger sizes of Elizabethan collar that can be the most problematic on a dog in the home. Problems such as the collar causing the dog to bang into doors and furniture, possibly injuring small children, or hitting the back of your legs, are significantly reduced. **Hound Surrounds™** keep the trusted cone shape, but the size of the cone has been reduced by 25% by overlapping and strongly riveting sections of the plastic. This allows the use of a proper size **Hound Surround™** while bringing the width down to what you would expect from a collar at least one size smaller. Also included are ventilation slots at the level of the ears to aid in the prevention and management of ear and skin infections.

Hound Surrounds™ come ready assembled with pre-constructed flat tabs for the dog's collar, making life easier for practice staff and allowing the **Hound Surround™** to lie comfortably on the dog's neck. Only three locking tabs are needed to size and fit a **Hound Surround™** to the dog and these tabs have been designed to be easy to use, but also flexible, strong and self locking. **Hound Surrounds™** are made from clear material, so maintaining the peripheral vision of the dog and come in three sizes, which correspond to conventional Elizabethan collar sizes: 25cm & 35cm. They also come in a Multipack of two of each size.



Harvey – with his size 30 Hound Surround

“Harvey with his “hoodie” on. This was marvellous, as you can see he took to this much better than the original one we were given for him.

There was much much less banging into things and people and he accepted it very well. Putting this on Harvey was easy as he accepted it without any fuss which in turn led to less anxiety on his behalf and also on ours.”

